

Saliva Ejectors



Always in a good
shape

If quality makes the difference



Leptrex Saliva Ejectors

have an excellent quality and provide an optimal suction.

They are extremely flexible and easily adaptable to any anatomical situation in the patient's mouth.

Once bent and placed in the mouth, the ejector maintains its shape due to the integrated solid metal wire, which also prevents breaching without kinking or even clogging.

A smooth tip allows air to flow even in the event of a blockage.

The tip can be removed if necessary. This allows even small particles to be aspirated from the oral cavity.



Relaxed treatment

Saliva Ejectors are used for almost every procedure. The working area is kept free from saliva, blood and other fluids and remains as dry as possible.



Position of the saliva ejector buccally




Position of the saliva ejector lingually

Two sizes, seven colours

Leptrex Saliva Ejectors are available in 13 cm and 15 cm length. So the perfect size for any situation can be chosen.

Choose the colour that matches to your practice best!



Images with friendly permission of Dr. M. Dobbertin, Frankfurt

How to order

+44 20 804 04912
info@leptrexltd.co.uk



Saliva Ejectors

Length: 13 cm
Content: 100 pcs.

Clear	Blue	Green	Yellow	Lavender	Orange	Pink
504-C	505-B	506-G	508-Y	510-L	511-O	512-P

Length: 15 cm
Content: 100 pcs.

Clear	Blue	Green	Yellow	Lavender	Orange	Pink
704-C	705-B	706-G	708-Y	710-L	711-O	712-P

The illustrations may show deviations from the original product in terms of colour etc.
Our general terms and conditions apply.



www.leptrexltd.co.uk